



SURGICRYL®

910


POLYGLACTINE



Type	Braided and coated multifilament
Composition	Polyglactine 910, a copolymer made of 90% glycolide and 10% L-lactide
Coating	Poly-glycolide-co-L-lactide and calcium stearate (<1%)
Colour	Violet
Tissue reaction	Minimal
Absorption	The hydrolytic action by which the material is broken down results in total absorption between 56 and 70 days. Approximately 50% of tensile strength remain after 21 days.
Presentation	Box with 12 sutures



*For quality
&
safety*



SURGICRYL®

910

POLYGLACTINE

CHARACTERISTICS

- High tensile strength
- Synthetic
- Breaks down by hydrolysis
- Braided structure

- Coated multifilament

- Hermetically sealed packing

BENEFITS

- Excellent tissue support
- Minimal tissue reaction
- Predictable absorption rate
- Excellent handling properties
 - High flexibility
 - High tensile strength
 - Excellent knot security
- Soft passage through the tissues with minimal sawing, tissue drag and trauma
- Good knot tie down / adjustability
- Reduced capillary action

- Guaranteed seal and product sterility

INDICATIONS General soft tissue closing and/or ligation especially in


- General surgery
- Skin closure
- Gastrointestinal surgery
- Plastic surgery
- Ophthalmic surgery
- Gynaecology
- Obstetrics
- Orthopaedics
- Urology

CONTRAINDICATIONS Cardiovascular surgery
Neurosurgery


Being absorbable not to use when an extended wound support is required.

Caution in patients with delayed wound healing (e.g. compromised immune system).

TENSILE STRENGTH



ABSORPTION



 Without needle

Metric	USP	Length		Order number
1.5	4/0	150 cm	<i>violet</i>	15150015
2	3/0	150 cm	<i>violet</i>	15200015
3	2/0	150 cm	<i>violet</i>	15300015
3.5	0	150 cm	<i>violet</i>	15350015
4	1	150 cm	<i>violet</i>	15400015
5	2	150 cm	<i>violet</i>	15500015
3	2/0	3 x 45 cm	<i>violet</i>	15300045-3
3.5	0	3 x 45 cm	<i>violet</i>	15350045-3
1.5	4/0	6 x 45 cm	<i>violet</i>	15150045-6
2	3/0	6 x 45 cm	<i>violet</i>	15200045-6
3	2/0	6 x 45 cm	<i>violet</i>	15300045-6
3.5	0	6 x 45 cm	<i>violet</i>	15350045-6
1.5	4/0	12 x 45 cm	<i>violet</i>	15150045
2	3/0	12 x 45 cm	<i>violet</i>	15200045
3	2/0	12 x 45 cm	<i>violet</i>	15300045
3.5	0	12 x 45 cm	<i>violet</i>	15350045
4	1	12 x 45 cm	<i>violet</i>	15400045
5	2	12 x 45 cm	<i>violet</i>	15500045
2	3/0	2 x 70 cm	<i>violet</i>	15200070-2
3	2/0	2 x 70 cm	<i>violet</i>	15300070-2
3.5	0	2 x 70 cm	<i>violet</i>	15350070-2
4	1	2 x 70 cm	<i>violet</i>	15400070-2
5	2	2 x 70 cm	<i>violet</i>	15500070-2
4	1	7 x 75 cm	<i>violet</i>	15400075

 Without needle

Ligapack



Metric	USP	Length		Order number
2	3/0	250 cm	<i>violet</i>	15200025
3	2/0	250 cm	<i>violet</i>	15300025
3.5	0	250 cm	<i>violet</i>	15350025
4	1	250 cm	<i>violet</i>	15400025
5	2	250 cm	<i>violet</i>	15500025



Round bodied taper point








1/2 circle



		Metric	USP	Length		Order number
HR 13 mm		0.7	6/0	75 cm	<i>violet</i>	15070113
		1	5/0	75 cm	<i>violet</i>	15100113
		1.5	4/0	75 cm	<i>violet</i>	15150113
HR 17 mm		0.7	6/0	75 cm	<i>violet</i>	15070117
		1	5/0	75 cm	<i>violet</i>	15100117
		1.5	4/0	75 cm	<i>violet</i>	15150117
		2	3/0	75 cm	<i>violet</i>	15200117
HR 22 mm		3	2/0	75 cm	<i>violet</i>	15300117
		1.5	4/0	75 cm	<i>violet</i>	15150122
		2	3/0	75 cm	<i>violet</i>	15200122
HR 26 mm		3	2/0	75 cm	<i>violet</i>	15300122
		4	1	75 cm	<i>violet</i>	15400122
		1	5/0	75 cm	<i>violet</i>	15100126
		1.5	4/0	75 cm	<i>violet</i>	15150126
HR 26 mm		2	3/0	75 cm	<i>violet</i>	15200126
		3	2/0	75 cm	<i>violet</i>	15300126
		3	2/0	75 cm	<i>undyed</i>	15300126U
		3.5	0	75 cm	<i>violet</i>	15350126
		4	1	75 cm	<i>violet</i>	15400126
		5	2	75 cm	<i>violet</i>	15500126
HR 30 mm		2	3/0	75 cm	<i>violet</i>	15200130
		3	2/0	75 cm	<i>violet</i>	15300130
		3.5	0	75 cm	<i>violet</i>	15350130
		4	1	75 cm	<i>violet</i>	15400130
HR 30 mm		5	2	90 cm	<i>violet</i>	15500130
		2	3/0	75 cm	<i>violet</i>	15200136
		3	2/0	75 cm	<i>violet</i>	15300136
HR 36 mm		3	2/0	75 cm	<i>undyed</i>	15300136U
		3.5	0	90 cm	<i>violet</i>	15350136
		4	1	90 cm	<i>violet</i>	15400136
		5	2	90 cm	<i>violet</i>	15500136








Round bodied taper point




		Metric	USP	Length		Order number
	HR 40 mm 	2	3/0	75 cm	<i>violet</i>	15200140
		3	2/0	75 cm	<i>violet</i>	15300140
		3.5	0	75 cm	<i>violet</i>	15350140-75
		3.5	0	90 cm	<i>violet</i>	15350140
		4	1	75 cm	<i>violet</i>	15400140-75
		4	1	90 cm	<i>violet</i>	15400140
		5	2	75 cm	<i>violet</i>	15500140-75
		5	2	90 cm	<i>violet</i>	15500140
	HR 45 mm 	4	1	90 cm	<i>violet</i>	15400145
		5	2	90 cm	<i>violet</i>	15500145
	HR 48 mm 	3	2/0	90 cm	<i>violet</i>	15300148
		3.5	0	90 cm	<i>violet</i>	15350148
		4	1	90 cm	<i>violet</i>	15400148
	HR 60 mm 	5	2	90 cm	<i>violet</i>	15500148
4		1	75 cm	<i>violet</i>	15400160	
HR 65 mm 	3	2/0	90 cm	<i>violet</i>	15300165	
	3.5	0	90 cm	<i>violet</i>	15350165	
HR 80 mm 	4	1	100 cm	<i>violet</i>	15400180	

















Round bodied taper point double armed


		Metric	USP	Length		Order number
	2 x HR 26 mm 	2	3/0	75 cm	<i>violet</i>	15200126-2
		3	2/0	75 cm	<i>violet</i>	15300126-2


 Round bodied taper point loop	Metric	USP	Length	Order number			
1/2 circle 	HR 45 mm		3.5	0	150 cm	<i>violet</i>	15350945












 Round bodied taper point	Metric	USP	Length	Order number			
5/8 circle 	FR 27 mm		3	2/0	75 cm	<i>violet</i>	15300327
			3.5	0	75 cm	<i>violet</i>	15350327







 Round bodied taper point	Metric	USP	Length	Order number			
3/8 circle 	DR 17 mm		1.5	4/0	45 cm	<i>violet</i>	15150217
	DR 22 mm		2	3/0	75 cm	<i>violet</i>	15200222
	DR 30 mm		4	1	75 cm	<i>violet</i>	15400230
	DR 32 mm		3	2/0	90 cm	<i>violet</i>	15300232
			3.5	0	75 cm	<i>violet</i>	15350232






 Reverse cutting	Metric	USP	Length	Order number			
1/2 circle 	HS 12 mm		0.7	6/0	45 cm	<i>violet</i>	15071412
	HS 17 mm		1	5/0	75 cm	<i>violet</i>	15101417
			1.5	4/0	75 cm	<i>violet</i>	15151417
	HS 22 mm		2	3/0	75 cm	<i>violet</i>	15201422
	HS 23 mm		5	2	90 cm	<i>violet</i>	15501423
			2	3/0	75 cm	<i>violet</i>	15201426
			3	2/0	75 cm	<i>violet</i>	15301426
			3.5	0	75 cm	<i>violet</i>	15351426
	HS 26 mm		4	1	75 cm	<i>violet</i>	15401426
			2	3/0	75 cm	<i>violet</i>	15201430
HS 30 mm							




Reverse cutting		Metric	USP	Length		Order number
1/2 circle 	HS 36 mm	3	2/0	75 cm	<i>violet</i>	15301436
		3.5	0	90 cm	<i>violet</i>	15351436
		4	1	90 cm	<i>violet</i>	15401436
		5	2	90 cm	<i>violet</i>	15501436
	HS 40 mm	2	3/0	90 cm	<i>violet</i>	15201440
		3	2/0	90 cm	<i>violet</i>	15301440
		3.5	0	90 cm	<i>violet</i>	15351440
		4	1	90 cm	<i>violet</i>	15401440
		5	2	90 cm	<i>violet</i>	15501440
	HS 48 mm	3.5	0	90 cm	<i>violet</i>	15351448
		4	1	90 cm	<i>violet</i>	15401448
		5	2	90 cm	<i>violet</i>	15501448




Reverse cutting		Metric	USP	Length		Order number
3/8 circle 	DS 12 mm	0.7	6/0	45 cm	<i>violet</i>	15071512
		1	5/0	45 cm	<i>violet</i>	15101512
		1.5	4/0	45 cm	<i>violet</i>	15151512
	DS 16 mm	0.7	6/0	45 cm	<i>violet</i>	15071516
		1	5/0	75 cm	<i>violet</i>	15101516
		1.5	4/0	75 cm	<i>violet</i>	15151516
		2	3/0	45 cm	<i>violet</i>	15201516-45
	DS 16 mm Easy Pass	2	3/0	50 cm	<i>undyed</i>	15201516UP
	DS 16 mm black	0.7	6/0	45 cm	<i>violet</i>	15071516VB
	DS 19 mm	1	5/0	75 cm	<i>violet</i>	15101519
		1.5	4/0	75 cm	<i>violet</i>	15151519
		2	3/0	75 cm	<i>violet</i>	15201519
	DS 19 mm Easy Pass	1.5	4/0	50 cm	<i>undyed</i>	15151519UP
		2	3/0	50 cm	<i>undyed</i>	15201519UP
	DS 19 mm black	1	5/0	75 cm	<i>violet</i>	15101519VB
		1.5	4/0	75 cm	<i>violet</i>	15151519VB
	DS 24 mm	1.5	4/0	75 cm	<i>violet</i>	15151524
		2	3/0	75 cm	<i>violet</i>	15201524
		3	2/0	75 cm	<i>violet</i>	15301524
		3.5	0	75 cm	<i>violet</i>	15351524




 Reverse cutting		Metric	USP	Length		Order number
	DS 24 mm Easy Pass 	2	3/0	50 cm	<i>undyed</i>	15201524UP
	DS 26 mm 	2	3/0	75 cm	<i>violet</i>	15201526
		3	2/0	75 cm	<i>violet</i>	15301526
	DS 26 mm Easy Pass 	2	3/0	75 cm	<i>undyed</i>	15201526UP
	DS 30 mm 	2	3/0	75 cm	<i>violet</i>	15201530
		3	2/0	75 cm	<i>violet</i>	15301530
		3	2/0	90 cm	<i>violet</i>	15301530-90
		3.5	0	75 cm	<i>violet</i>	15351530
		3.5	0	90 cm	<i>violet</i>	15351530-90
	DS 30 mm Easy Pass 	3	2/0	75 cm	<i>undyed</i>	15301530UP-75
	DS 36 mm 	3	2/0	75 cm	<i>violet</i>	15301536
		3.5	0	75 cm	<i>violet</i>	15351536
	DS 40 mm 	4	1	75 cm	<i>violet</i>	15401536
		3.5	0	75 cm	<i>violet</i>	15351540
	DS 45 mm 	4	1	75 cm	<i>violet</i>	15401540
3		2/0	75 cm	<i>violet</i>	15301545	
DS 80 mm 	7	5	150 cm	<i>violet</i>	15701580	




 Round bodied taper cutting point		Metric	USP	Length		Order number
	HRT 13 mm 	1.5	4/0	75 cm	<i>violet</i>	15150513
	HRT 17 mm 	1.5	4/0	75 cm	<i>violet</i>	15150517
		2	3/0	75 cm	<i>violet</i>	15200517
	HRT 26 mm 	1.5	4/0	75 cm	<i>violet</i>	15150526
		2	3/0	75 cm	<i>violet</i>	15200526
		3	2/0	75 cm	<i>violet</i>	15300526
	HRT 30 mm 	3	2/0	75 cm	<i>violet</i>	15300530
		3.5	0	75 cm	<i>violet</i>	15350530
	4	1	75 cm	<i>violet</i>	15400530	






 Round bodied taper cutting point	Metric	USP	Length	Order number	
 1/2 circle	 HRT 36 mm	3	2/0	90 cm	<i>violet</i> 15300536
		3.5	0	90 cm	<i>violet</i> 15350536
		4	1	90 cm	<i>violet</i> 15400536
	 HRT 40 mm	3.5	0	90 cm	<i>violet</i> 15350540
		4	1	90 cm	<i>violet</i> 15400540
		5	2	90 cm	<i>violet</i> 15500540
	 HRT 48 mm	8	6	150 cm	<i>violet</i> 15800548





 Round bodied taper cutting point	Metric	USP	Length	Order number	
 3/8 circle	 DRT 17 mm	1.5	4/0	75 cm	<i>violet</i> 15150417
		2	3/0	75 cm	<i>violet</i> 15200417




 Round bodied taper cutting point	Metric	USP	Length	Order number
 3/8 circle	 DRT 13 mm	1	75 cm	<i>violet</i> 15100530-0413








 Round bodied taper point	Metric	USP	Length	Order number
 1/2 circle	 HR 40 mm	4	180 cm	<i>violet</i> 15400140-0540-180








 Cutting	Metric	USP	Length	Order number	
 Straight	 GS 60 mm	2	3/0	75 cm	<i>violet</i> 15202160
		3	2/0	75 cm	<i>violet</i> 15302160





 Blunt		Metric	USP	Length		Order number
	HRB 36 mm	 3.5	0	90 cm	<i>violet</i>	15351836
		4	1	90 cm	<i>violet</i>	15401836
	HRB 40 mm	 4	1	90 cm	<i>violet</i>	15401840
	HRB 45 mm	 5	2	90 cm	<i>violet</i>	15501845








 Blunt		Metric	USP	Length		Order number
	DRB 63 mm	 4	1	100 cm	<i>violet</i>	15401963
	DRB 85 mm	 5	2	75 cm	<i>violet</i>	15501985

 Micro needle taper point		Metric	USP	Length		Order number
	MDR 6.4 mm 150µm	 0.4	8/0	20 cm	<i>violet</i>	M910087

 Side cutting spatula double armed	Metric	USP	Length	Order number			
1/2 circle 	2 x HSP 8 mm 200 µm		0.4	8/0	45 cm	<i>violet</i>	O910081
	2 x HSP 8 mm 200 µm		0.5	7/0	45 cm	<i>violet</i>	O91007
	2 x HSP 8 mm 300 µm		0.7	6/0	45 cm	<i>violet</i>	O910061
	2 x HSP 8 mm 300 µm		0.7	6/0	45 cm	<i>undyed</i>	O910061U
2 x HSP 8 mm 350 µm		1	5/0	45 cm	<i>violet</i>	O910052	

 Side cutting spatula double armed	Metric	USP	Length	Order number			
3/8 circle 	2 x DSP 6 mm 200 µm		0.4	8/0	45 cm	<i>violet</i>	O91008
	2 x DSP 6 mm 220 µm		0.5	7/0	45 cm	<i>violet</i>	O910072
	2 x DSP 6.5 mm 200 µm		0.4	8/0	30 cm	<i>violet</i>	O910082
	2 x DSP 6.5 mm 200 µm		0.5	7/0	30 cm	<i>violet</i>	O910077
	2 x DSP 11 mm 350 µm		1	5/0	45 cm	<i>violet</i>	O910051

 Side cutting spatula	Metric	USP	Length	Order number			
3/8 circle 	DSP 6.5 mm 200 µm		0.4	8/0	30 cm	<i>violet</i>	O910082-1
	DSP 6.5 mm 200 µm		0.5	7/0	30 cm	<i>violet</i>	O910077-1

 Side cutting spatula double armed	Metric	USP	Length	Order number			
1/4 circle 	2 x VSP 8 mm 300 µm		0.7	6/0	45 cm	<i>violet</i>	O91006
	2 x VSP 8 mm 300 µm		0.7	6/0	20 cm	<i>undyed</i>	O91006U-20
	2 x VSP 8 mm 300 µm		0.7	6/0	45 cm	<i>undyed</i>	O91006U
	2 x VSP 8 mm 350 µm		1	5/0	45 cm	<i>violet</i>	O91005
	2 x VSP 8 mm 350 µm		1	5/0	45 cm	<i>undyed</i>	O91005U