


# SURGICRYL<sup>®</sup>

910


POLYGLACTINE



<b>Type</b>	Braided and coated multifilament
<b>Composition</b>	Polyglactine 910, a copolymer made of 90% glycolide and 10% L-lactide
<b>Coating</b>	Poly-glycolide-co-L-lactide and calcium stearate (<1%)
<b>Colour</b>	Violet
<b>Tissue reaction</b>	Minimal
<b>Absorption</b>	The hydrolytic action by which the material is broken down results in total absorption between 56 and 70 days. Approximately 50% of tensile strength remain after 21 days.
<b>Presentation</b>	Box with 12 sutures; Cassette



*For quality  
&  
safety*



[www.sutures.be](http://www.sutures.be)

# SURGICRYL®

## 910

### POLYGLACTINE

### CHARACTERISTICS

- High tensile strength
- Synthetic
- Breaks down by hydrolysis
- Braided structure
  
- Coated multifilament
  
- Hermetically sealed packing

### BENEFITS

- Excellent tissue support
- Minimal tissue reaction
- Predictable absorption rate
- Excellent handling properties
  - High flexibility
  - High tensile strength
  - Excellent knot security
- Soft passage through the tissues with minimal sawing, tissue drag and trauma
- Good knot tie down / adjustability
- Reduced capillary action
  
- Guaranteed seal and product sterility

**INDICATIONS** General soft tissue closing and/or ligation especially in


- General surgery
- Skin closure
- Gastrointestinal surgery
- Plastic surgery
- Ophthalmic surgery
- Gynaecology
- Obstetrics
- Orthopaedics
- Urology

**CONTRAINDICATIONS** Cardiovascular surgery  
Neurosurgery


Being absorbable not to use when an extended wound support is required.

Caution in patients with delayed wound healing (e.g. compromised immune system).

### TENSILE STRENGTH



### ABSORPTION



○ Without needle

Metric	USP	Length		Order number
1.5	4/0	150 cm	<i>violet</i>	15150015
2	3/0	150 cm	<i>violet</i>	15200015
3	2/0	150 cm	<i>violet</i>	15300015
3.5	0	150 cm	<i>violet</i>	15350015
4	1	150 cm	<i>violet</i>	15400015
5	2	150 cm	<i>violet</i>	15500015
3	2/0	3 x 45 cm	<i>violet</i>	15300045-3
3.5	0	3 x 45 cm	<i>violet</i>	15350045-3
1.5	4/0	6 x 45 cm	<i>violet</i>	15150045-6
2	3/0	6 x 45 cm	<i>violet</i>	15200045-6
3	2/0	6 x 45 cm	<i>violet</i>	15300045-6
3.5	0	6 x 45 cm	<i>violet</i>	15350045-6
1.5	4/0	12 x 45 cm	<i>violet</i>	15150045
2	3/0	12 x 45 cm	<i>violet</i>	15200045
3	2/0	12 x 45 cm	<i>violet</i>	15300045
3.5	0	12 x 45 cm	<i>violet</i>	15350045
4	1	12 x 45 cm	<i>violet</i>	15400045
5	2	12 x 45 cm	<i>violet</i>	15500045
2	3/0	2 x 70 cm	<i>violet</i>	15200070-2
3	2/0	2 x 70 cm	<i>violet</i>	15300070-2
3.5	0	2 x 70 cm	<i>violet</i>	15350070-2
4	1	2 x 70 cm	<i>violet</i>	15400070-2
5	2	2 x 70 cm	<i>violet</i>	15500070-2
4	1	7 x 75 cm	<i>violet</i>	15400075

○ Without needle

Ligapack



Metric	USP	Length		Order number
2	3/0	250 cm	<i>violet</i>	15200025
3	2/0	250 cm	<i>violet</i>	15300025
3.5	0	250 cm	<i>violet</i>	15350025
4	1	250 cm	<i>violet</i>	15400025
5	2	250 cm	<i>violet</i>	15500025



## Round bodied taper point








1/2 circle



		Metric	USP	Length		Order number
HR 13 mm		0.7	6/0	75 cm	<i>violet</i>	15070113
		1	5/0	75 cm	<i>violet</i>	15100113
		1.5	4/0	75 cm	<i>violet</i>	15150113
HR 17 mm		0.7	6/0	75 cm	<i>violet</i>	15070117
		1	5/0	75 cm	<i>violet</i>	15100117
		1.5	4/0	75 cm	<i>violet</i>	15150117
		2	3/0	75 cm	<i>violet</i>	15200117
HR 22 mm		3	2/0	75 cm	<i>violet</i>	15300117
		1.5	4/0	75 cm	<i>violet</i>	15150122
		2	3/0	75 cm	<i>violet</i>	15200122
HR 26 mm		3	2/0	75 cm	<i>violet</i>	15300122
		4	1	75 cm	<i>violet</i>	15400122
		1	5/0	75 cm	<i>violet</i>	15100126
		1.5	4/0	75 cm	<i>violet</i>	15150126
HR 30 mm		2	3/0	75 cm	<i>violet</i>	15200126
		3	2/0	75 cm	<i>violet</i>	15300126
		3	2/0	75 cm	<i>undyed</i>	15300126U
		3.5	0	75 cm	<i>violet</i>	15350126
		4	1	75 cm	<i>violet</i>	15400126
		5	2	75 cm	<i>violet</i>	15500126
		2	3/0	75 cm	<i>violet</i>	15200130
HR 36 mm		3	2/0	75 cm	<i>violet</i>	15300130
		3.5	0	75 cm	<i>violet</i>	15350130
		4	1	75 cm	<i>violet</i>	15400130
		5	2	90 cm	<i>violet</i>	15500130
HR 36 mm		2	3/0	75 cm	<i>violet</i>	15200136
		3	2/0	75 cm	<i>violet</i>	15300136
		3	2/0	75 cm	<i>undyed</i>	15300136U
		3.5	0	90 cm	<i>violet</i>	15350136
		4	1	90 cm	<i>violet</i>	15400136
		5	2	90 cm	<i>violet</i>	15500136








### Round bodied taper point




		Metric	USP	Length		Order number
	HR 40 mm 	2	3/0	75 cm	<i>violet</i>	15200140
		3	2/0	75 cm	<i>violet</i>	15300140
		3.5	0	75 cm	<i>violet</i>	15350140-75
		3.5	0	90 cm	<i>violet</i>	15350140
		4	1	75 cm	<i>violet</i>	15400140-75
		4	1	90 cm	<i>violet</i>	15400140
		5	2	75 cm	<i>violet</i>	15500140-75
		5	2	90 cm	<i>violet</i>	15500140
	HR 45 mm 	4	1	90 cm	<i>violet</i>	15400145
		5	2	90 cm	<i>violet</i>	15500145
	HR 48 mm 	3	2/0	90 cm	<i>violet</i>	15300148
		3.5	0	90 cm	<i>violet</i>	15350148
		4	1	90 cm	<i>violet</i>	15400148
	HR 60 mm 	5	2	90 cm	<i>violet</i>	15500148
4		1	75 cm	<i>violet</i>	15400160	
HR 65 mm 	3	2/0	90 cm	<i>violet</i>	15300165	
	3.5	0	90 cm	<i>violet</i>	15350165	
HR 80 mm 	4	1	100 cm	<i>violet</i>	15400180	

















### Round bodied taper point double armed


		Metric	USP	Length		Order number
	2 x HR 26 mm 	2	3/0	75 cm	<i>violet</i>	15200126-2
		3	2/0	75 cm	<i>violet</i>	15300126-2


 Round bodied taper point loop	Metric	USP	Length	Order number			
<b>1/2 circle</b> 	HR 45 mm		3.5	0	150 cm	<i>violet</i>	15350945












 Round bodied taper point	Metric	USP	Length	Order number			
<b>5/8 circle</b> 	FR 27 mm		3	2/0	75 cm	<i>violet</i>	15300327
			3.5	0	75 cm	<i>violet</i>	15350327







 Round bodied taper point	Metric	USP	Length	Order number			
<b>3/8 circle</b> 	DR 17 mm		1.5	4/0	45 cm	<i>violet</i>	15150217
	DR 22 mm		2	3/0	75 cm	<i>violet</i>	15200222
	DR 30 mm		4	1	75 cm	<i>violet</i>	15400230
	DR 32 mm		3	2/0	90 cm	<i>violet</i>	15300232
			3.5	0	75 cm	<i>violet</i>	15350232

 Reverse cutting	Metric	USP	Length	Order number			
<b>1/2 circle</b> 	HS 12 mm		0.7	6/0	45 cm	<i>violet</i>	15071412
	HS 17 mm		1	5/0	75 cm	<i>violet</i>	15101417
			1.5	4/0	75 cm	<i>violet</i>	15151417
	HS 22 mm		2	3/0	75 cm	<i>violet</i>	15201422
	HS 23 mm		5	2	90 cm	<i>violet</i>	15501423
			2	3/0	75 cm	<i>violet</i>	15201426
			3	2/0	75 cm	<i>violet</i>	15301426
			3.5	0	75 cm	<i>violet</i>	15351426
	HS 26 mm		4	1	75 cm	<i>violet</i>	15401426
			2	3/0	75 cm	<i>violet</i>	15201430
HS 30 mm							






Reverse cutting		Metric	USP	Length		Order number
<b>1/2 circle</b> 	HS 36 mm	3	2/0	75 cm	<i>violet</i>	15301436
		3.5	0	90 cm	<i>violet</i>	15351436
		4	1	90 cm	<i>violet</i>	15401436
		5	2	90 cm	<i>violet</i>	15501436
	HS 40 mm	2	3/0	90 cm	<i>violet</i>	15201440
		3	2/0	90 cm	<i>violet</i>	15301440
		3.5	0	90 cm	<i>violet</i>	15351440
		4	1	90 cm	<i>violet</i>	15401440
		5	2	90 cm	<i>violet</i>	15501440
	HS 48 mm	3.5	0	90 cm	<i>violet</i>	15351448
		4	1	90 cm	<i>violet</i>	15401448
		5	2	90 cm	<i>violet</i>	15501448




Reverse cutting		Metric	USP	Length		Order number
<b>3/8 circle</b> 	DS 12 mm	0.7	6/0	45 cm	<i>violet</i>	15071512
		1	5/0	45 cm	<i>violet</i>	15101512
		1.5	4/0	45 cm	<i>violet</i>	15151512
	DS 16 mm	0.7	6/0	45 cm	<i>violet</i>	15071516
		1	5/0	75 cm	<i>violet</i>	15101516
		1.5	4/0	75 cm	<i>violet</i>	15151516
		2	3/0	45 cm	<i>violet</i>	15201516-45
	DS 16 mm Easy Pass	2	3/0	50 cm	<i>undyed</i>	15201516UP
	DS 16 mm black	0.7	6/0	45 cm	<i>violet</i>	15071516VB
	DS 19 mm	1	5/0	75 cm	<i>violet</i>	15101519
		1.5	4/0	75 cm	<i>violet</i>	15151519
		2	3/0	75 cm	<i>violet</i>	15201519
	DS 19 mm Easy Pass	1.5	4/0	50 cm	<i>undyed</i>	15151519UP
		2	3/0	50 cm	<i>undyed</i>	15201519UP
	DS 19 mm black	1	5/0	75 cm	<i>violet</i>	15101519VB
		1.5	4/0	75 cm	<i>violet</i>	15151519VB
	DS 24 mm	1.5	4/0	75 cm	<i>violet</i>	15151524
		2	3/0	75 cm	<i>violet</i>	15201524
3		2/0	75 cm	<i>violet</i>	15301524	
	3.5	0	75 cm	<i>violet</i>	15351524	




 Reverse cutting		Metric	USP	Length		Order number
	DS 24 mm Easy Pass 	2	3/0	50 cm	<i>undyed</i>	15201524UP
	DS 26 mm 	2	3/0	75 cm	<i>violet</i>	15201526
		3	2/0	75 cm	<i>violet</i>	15301526
	DS 26 mm Easy Pass 	2	3/0	75 cm	<i>undyed</i>	15201526UP
	DS 30 mm 	2	3/0	75 cm	<i>violet</i>	15201530
		3	2/0	75 cm	<i>violet</i>	15301530
		3	2/0	90 cm	<i>violet</i>	15301530-90
		3.5	0	75 cm	<i>violet</i>	15351530
		3.5	0	90 cm	<i>violet</i>	15351530-90
	DS 30 mm Easy Pass 	3	2/0	75 cm	<i>undyed</i>	15301530UP-75
	DS 36 mm 	3	2/0	75 cm	<i>violet</i>	15301536
		3.5	0	75 cm	<i>violet</i>	15351536
		4	1	75 cm	<i>violet</i>	15401536
	DS 40 mm 	3.5	0	75 cm	<i>violet</i>	15351540
		4	1	75 cm	<i>violet</i>	15401540
DS 45 mm 	3	2/0	75 cm	<i>violet</i>	15301545	
DS 80 mm 	7	5	150 cm	<i>violet</i>	15701580	





 Round bodied taper cutting point		Metric	USP	Length		Order number
	HRT 13 mm 	1.5	4/0	75 cm	<i>violet</i>	15150513
	HRT 17 mm 	1.5	4/0	75 cm	<i>violet</i>	15150517
		2	3/0	75 cm	<i>violet</i>	15200517
	HRT 26 mm 	1.5	4/0	75 cm	<i>violet</i>	15150526
		2	3/0	75 cm	<i>violet</i>	15200526
		3	2/0	75 cm	<i>violet</i>	15300526
		3	2/0	75 cm	<i>violet</i>	15300530
	HRT 30 mm 	3.5	0	75 cm	<i>violet</i>	15350530
		4	1	75 cm	<i>violet</i>	15400530











 Round bodied taper cutting point	Metric	USP	Length	Order number	
 <b>1/2 circle</b>	HRT 36 mm 	3	2/0	90 cm	<i>violet</i> 15300536
		3.5	0	90 cm	<i>violet</i> 15350536
		4	1	90 cm	<i>violet</i> 15400536
	HRT 40 mm 	3.5	0	90 cm	<i>violet</i> 15350540
		4	1	90 cm	<i>violet</i> 15400540
		5	2	90 cm	<i>violet</i> 15500540
	HRT 48 mm 	8	6	150 cm	<i>violet</i> 15800548





 Round bodied taper cutting point	Metric	USP	Length	Order number	
 <b>3/8 circle</b>	DRT 17 mm 	1.5	4/0	75 cm	<i>violet</i> 15150417
		2	3/0	75 cm	<i>violet</i> 15200417




 Round bodied taper cutting point	Metric	USP	Length	Order number
 <b>3/8 circle</b>	DRT 13 mm 	1	75 cm	<i>violet</i> 15100530-0413








 Round bodied taper point	Metric	USP	Length	Order number
 <b>1/2 circle</b>	HR 40 mm 	4	180 cm	<i>violet</i> 15400140-0540-180
	HRT 40 mm 			








 Cutting	Metric	USP	Length	Order number
 <b>Straight</b>	GS 60 mm 	2	75 cm	<i>violet</i> 15202160
		3	75 cm	<i>violet</i> 15302160





 Blunt		Metric	USP	Length		Order number
	HRB 36 mm	 3.5	0	90 cm	<i>violet</i>	15351836
		4	1	90 cm	<i>violet</i>	15401836
	HRB 40 mm	 4	1	90 cm	<i>violet</i>	15401840
	HRB 45 mm	 5	2	90 cm	<i>violet</i>	15501845








 Blunt		Metric	USP	Length		Order number
	DRB 63 mm	 4	1	100 cm	<i>violet</i>	15401963
	DRB 85 mm	 5	2	75 cm	<i>violet</i>	15501985

 Micro needle taper point		Metric	USP	Length		Order number
	MDR 6.4 mm 150µm	 0.4	8/0	20 cm	<i>violet</i>	M910087

 Side cutting spatula double armed	Metric	USP	Length	Order number			
<b>1/2 circle</b> 	2 x HSP 8 mm 200 µm		0.4	8/0	45 cm	<i>violet</i>	O910081
	2 x HSP 8 mm 200 µm		0.5	7/0	45 cm	<i>violet</i>	O91007
	2 x HSP 8 mm 300 µm		0.7	6/0	45 cm	<i>violet</i>	O910061
	2 x HSP 8 mm 300 µm		0.7	6/0	45 cm	<i>undyed</i>	O910061U
2 x HSP 8 mm 350 µm		1	5/0	45 cm	<i>violet</i>	O910052	

 Side cutting spatula double armed	Metric	USP	Length	Order number			
<b>3/8 circle</b> 	2 x DSP 6 mm 200 µm		0.4	8/0	45 cm	<i>violet</i>	O91008
	2 x DSP 6 mm 220 µm		0.5	7/0	45 cm	<i>violet</i>	O910072
	2 x DSP 6.5 mm 200 µm		0.4	8/0	30 cm	<i>violet</i>	O910082
	2 x DSP 6.5 mm 200 µm		0.5	7/0	30 cm	<i>violet</i>	O910077
	2 x DSP 11 mm 350 µm		1	5/0	45 cm	<i>violet</i>	O910051

 Side cutting spatula	Metric	USP	Length	Order number			
<b>3/8 circle</b> 	DSP 6.5 mm 200 µm		0.4	8/0	30 cm	<i>violet</i>	O910082-1
	DSP 6.5 mm 200 µm		0.5	7/0	30 cm	<i>violet</i>	O910077-1

 Side cutting spatula double armed	Metric	USP	Length	Order number			
<b>1/4 circle</b> 	2 x VSP 8 mm 300 µm		0.7	6/0	45 cm	<i>violet</i>	O91006
	2 x VSP 8 mm 300 µm		0.7	6/0	20 cm	<i>undyed</i>	O91006U-20
	2 x VSP 8 mm 300 µm		0.7	6/0	45 cm	<i>undyed</i>	O91006U
	2 x VSP 8 mm 350 µm		1	5/0	45 cm	<i>violet</i>	O91005
	2 x VSP 8 mm 350 µm		1	5/0	45 cm	<i>undyed</i>	O91005U

## IN CASSETTES

Use within 4 months after first opening

Metric	USP	Length	Order number
1.5	4/0	15 m	17151
2	3/0	15 m	17201
2	3/0	25 m	17202
3	2/0	15 m	17301
3	2/0	25 m	17302
3	2/0	50 m	17305
3.5	0	15 m	17351
3.5	0	25 m	17352
3.5	0	50 m	17355
4	1	15 m	17401
4	1	25 m	17402
4	1	50 m	17405
5	2	15 m	17501
5	2	25 m	17502
5	2	50 m	17505
5	2	100 m	175010
6	3&4	15 m	17601
6	3&4	25 m	17602
6	3&4	50 m	17605
7	5	15 m	17701
7	5	25 m	17702
7	5	50 m	17705
8	6	15 m	17801
8	6	25 m	17802
8	6	50 m	17805

