



SURGICRYL[®]

910


POLYGLACTINE



Type	Braided and coated multifilament
Composition	Polyglactine 910, a copolymer made of 90% glycolide and 10% L-lactide
Coating	Poly-glycolide-co-L-lactide and calcium stearate (<1%)
Colour	Violet
Tissue reaction	Minimal
Absorption	The hydrolytic action by which the material is broken down results in total absorption between 56 and 70 days. Approximately 50% of tensile strength remain after 21 days.
Presentation	Box with 12 sutures



*For quality
&
safety*



www.sutures.be

SURGICRYL®

910

POLYGLACTINE

CHARACTERISTICS

- High tensile strength
- Synthetic
- Breaks down by hydrolysis
- Braided structure

- Coated multifilament

- Hermetically sealed packing

BENEFITS

- Excellent tissue support
- Minimal tissue reaction
- Predictable absorption rate
- Excellent handling properties
 - High flexibility
 - High tensile strength
 - Excellent knot security
- Soft passage through the tissues with minimal sawing, tissue drag and trauma
- Good knot tie down / adjustability
- Reduced capillary action

- Guaranteed seal and product sterility

INDICATIONS General soft tissue closing and/or ligation especially in


- General surgery
- Skin closure
- Gastrointestinal surgery
- Plastic surgery
- Ophthalmic surgery
- Gynaecology
- Obstetrics
- Orthopaedics
- Urology

CONTRAINDICATIONS Cardiovascular surgery
Neurosurgery


Being absorbable not to use when an extended wound support is required.






Caution in patients with delayed wound healing (e.g. compromised immune system).





TENSILE STRENGTH












ABSORPTION







 Round bodied taper point		Metric	USP	Length		Order number
	HR 13 mm 	0.7	6/0	75 cm	<i>violet</i>	15070113
		1	5/0	75 cm	<i>violet</i>	15100113
		1.5	4/0	75 cm	<i>violet</i>	15150113
	HR 17 mm 	0.7	6/0	75 cm	<i>violet</i>	15070117
		1	5/0	75 cm	<i>violet</i>	15100117
		1.5	4/0	75 cm	<i>violet</i>	15150117
		2	3/0	75 cm	<i>violet</i>	15200117
		3	2/0	75 cm	<i>violet</i>	15300117
	HR 22 mm 	1.5	4/0	75 cm	<i>violet</i>	15150122
2		3/0	75 cm	<i>violet</i>	15200122	
3		2/0	75 cm	<i>violet</i>	15300122	

 Round bodied taper point		Metric	USP	Length		Order number
	DR 17 mm 	1.5	4/0	45 cm	<i>violet</i>	15150217
	DR 22 mm 	2	3/0	75 cm	<i>violet</i>	15200222

 Reverse cutting		Metric	USP	Length		Order number
	HS 12 mm 	0.7	6/0	45 cm	<i>violet</i>	15071412
	HS 17 mm 	1	5/0	75 cm	<i>violet</i>	15101417
		1.5	4/0	75 cm	<i>violet</i>	15151417
	HS 22 mm 	2	3/0	75 cm	<i>violet</i>	15201422
	HS 23 mm 	5	2	90 cm	<i>violet</i>	15501423
	HS 26 mm 	2	3/0	75 cm	<i>violet</i>	15201426
		3	2/0	75 cm	<i>violet</i>	15301426

 Reverse cutting		Metric	USP	Length		Order number
	DS 12 mm	0.7	6/0	45 cm	<i>violet</i>	15071512
		1	5/0	45 cm	<i>violet</i>	15101512
		1.5	4/0	45 cm	<i>violet</i>	15151512
	DS 16 mm	0.7	6/0	45 cm	<i>violet</i>	15071516
		1	5/0	75 cm	<i>violet</i>	15101516
		1.5	4/0	75 cm	<i>violet</i>	15151516
		2	3/0	45 cm	<i>violet</i>	15201516-45
	DS 16 mm Easy Pass	2	3/0	50 cm	<i>undyed</i>	15201516UP
	DS 16 mm black	0.7	6/0	45 cm	<i>violet</i>	15071516VB
	DS 19 mm	1	5/0	75 cm	<i>violet</i>	15101519
		1.5	4/0	75 cm	<i>violet</i>	15151519
		2	3/0	75 cm	<i>violet</i>	15201519
	DS 19 mm Easy Pass	1.5	4/0	50 cm	<i>undyed</i>	15151519UP
		2	3/0	50 cm	<i>undyed</i>	15201519UP
	DS 19 mm black	1	5/0	75 cm	<i>violet</i>	15101519VB
		1.5	4/0	75 cm	<i>violet</i>	15151519VB
	DS 24 mm	1.5	4/0	75 cm	<i>violet</i>	15151524
		2	3/0	75 cm	<i>violet</i>	15201524
3		2/0	75 cm	<i>violet</i>	15301524	
DS 24 mm Easy Pass	2	3/0	50 cm	<i>undyed</i>	15201524UP	

 Round bodied taper cutting point		Metric	USP	Length		Order number
	HRT 13 mm	1.5	4/0	75 cm	<i>violet</i>	15150513
	HRT 17 mm	1.5	4/0	75 cm	<i>violet</i>	15150517
		2	3/0	75 cm	<i>violet</i>	15200517

 Round bodied taper cutting point		Metric	USP	Length		Order number
	DRT 17 mm	1.5	4/0	75 cm	<i>violet</i>	15150417
		2	3/0	75 cm	<i>violet</i>	15200417